ENGLISH:**I. INTENDED USE:**

The Fentanyl Analogues 2 ng/mL, whole blood toxicology control is for use as a control material to evaluate the performance of an analytical test method.

II. SUMMARY AND PRINCIPLES:

Several different techniques are used for evaluating or estimating the variance of results. The three subjects summarized below must be considered with any test method.

1. PREVENTIVE MEASURES:

These measures are usually contained in the design of the test method and include consideration for reagents, equipment, and operator errors. These measures are designed to minimize variance.

2. CONTROL MEASURES:

When a control sample is analyzed at the same time and in the same manner as a patient specimen, an estimate of variance is obtained for the test method. This estimate of variance can be compared to the acceptable limits of variance of the test method.

3. STATISTICAL ANALYSIS OF PATIENT RESULTS:

As an aid in evaluating overall test results, the past experience of expected results can be compared to the results of any given test run. For example, it would not be expected that all results of a given test run be in an elevated range.

Control materials are widely used as a means to aid in the evaluation of test results. The following subjects are to be considered in the use of any control material.

1. Multi-Level
2. Matrix
3. Availability
4. Form
5. Variety

NORMAL / ELEVATED
HUMAN / ANIMAL / CHEMICAL
SUFFICIENT FOR STATISTICS
LIQUID / FROZEN / DRIED
DIFFERENT THAN CALIBRATORS

III. PRECAUTIONS:

1. Although the whole blood has been tested and found negative for HBsAg by RIA and HIV by EIA, the control material should be treated as any other potentially infectious agent.
2. For *in vitro* diagnostic use only.

IV. STORAGE AND STABILITY:

1. Store fresh frozen control material at or below -10°C (14°F). Stable to expiration date printed on the insert and label.
2. Store thawed control material at or below 2-8°C (35-46°F). See table for thawed stability.

V. PROCEDURE:

1. Allow fresh frozen control material to thaw at room temperature with cap on.
2. Swirl gently 3-4 minutes to ensure a homogeneous mixture.
3. Swirl gently each time an aliquot is removed to ensure a homogeneous mixture.
4. Assay control material in same manner as patient specimens, following the exact same instructions from the entire test method.
5. Record the results obtained on a control chart that describes statistical limits for the test method and the particular lot of control material.

VI. LIMITATIONS:

1. Results are dependent upon proper storage and adequate mixing.
2. Control material approximates a patient specimen.

VII. EXPECTED VALUES:

1. Laboratories should establish their own mean values; an individual laboratory's mean of several determinations may not duplicate the values listed below.

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REF 10189

Product Description: Fentanyl Analogues 2 ng/mL

Product Number:

LOT C7598

Lot Number:

 2022-11-30

Expiration Date:

Analyte	Stability(Days)	Gravimetrically Weighed Value	Verified Value	Units
4-ANPP	*N/A	2	1.7	ng/mL
Acetylfentanyl	*N/A	2	1.5	ng/mL
Acrylfentanyl	*N/A	2	1.7	ng/mL
Alfentanil	25	2	2.1	ng/mL
Butyryl Fentanyl	25	2	2.1	ng/mL
Carfentanyl	*N/A	2	1.9	ng/mL
cis-3-Methylfentanyl	*N/A	2	1.9	ng/mL
Crotonylfentanyl	25	2	*Present	ng/mL
Fentanyl	25	2	1.7	ng/mL
Furanyl Fentanyl	25	2	2.0	ng/mL
Isobutyrylfentanyl	25	2	1.9	ng/mL
m-Methylmethoxyacetylfentanyl	25	2	2.0	ng/mL
Methoxyacetylfentanyl	25	2	2.1	ng/mL
Norfentanyl	25	2	2.2	ng/mL
Ocfentanyl	25	2	*Present	ng/mL
ortho-Fluorofentanyl	*N/A	2	1.7	ng/mL
para-Fluorobutyrylfentanyl	*N/A	2	1.8	ng/mL
para-Fluoroisobutyrylfentanyl	*N/A	2	1.9	ng/mL
para-Fluorofentanyl	*N/A	2	2.0	ng/mL
Sufentanil	25	2	1.4	ng/mL
Tetrahydrofuran Fentanyl	25	2	2.3	ng/mL
U-47700	25	2	1.7	ng/mL
U-49900	25	2	1.7	ng/mL
U-51754	25	2	2.0	ng/mL
Valeryl fentanyl	25	2	1.8	ng/mL
Cyclopropylfentanyl	25	2	2.1	ng/mL

Matrix: Whole Blood

Product Form: Frozen

Storage: -10°C (14°F)

Vial Volume: 3 mL

*Analytes have been weighed and spiked into the control; however, no claim will be made for expected values or stability.

UTAK's express and implied warranties (including merchantability and fitness) are conditioned on the observance of UTAK's insert directions with respect to the use of UTAK's products.

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For technical support please call: UTAK Technical Service (888) 882-5522

Page | 4 of 5

IVD

UTAK Laboratories, Inc.
25020 Avenue Tibbitts
Valencia, CA 91355
T: (888) 882-5522
F: (661) 294-9272
E: welovecontrol@utak.com

EMERGO EUROPE
PRINSESSEGRACHT 20
2514 AP THE HAGUE
THE NETHERLANDS