ENGLISH:

I. INTENDED USE:

The Benzodiazepines Plus 100 ng/mL, urine toxicology control is for use as a control material to evaluate the performance of an analytical test method.

II. SUMMARY AND PRINCIPLES:

Several different techniques are used for evaluating or estimating the variance of results. The three subjects summarized below must be considered with any test method.

1. PREVENTIVE MEASURES:

These measures are usually contained in the design of the test method and include consideration for reagents, equipment, and operator errors. These measures are designed to minimize variance.

2. CONTROL MEASURES:

When a control sample is analyzed at the same time and in the same manner as a patient specimen, an estimate of variance is obtained for the test method. This estimate of variance can be compared to the acceptable limits of variance of the test method.

3. STATISTICAL ANALYSIS OF PATIENT RESULTS:

As an aid in evaluating overall test results, the past experience of expected results can be compared to the results of any given test run. For example, it would not be expected that all results of a given test run be in an elevated range.

Control materials are widely used as a means to aid in the evaluation of test results. The following subjects are to be considered in the use of any control material.

- | | |
|-----------------|----------------------------|
| 1. Multi-Level | NORMAL / ELEVATED |
| 2. Matrix | HUMAN / ANIMAL / CHEMICAL |
| 3. Availability | SUFFICIENT FOR STATISTICS |
| 4. Form | LIQUID / FROZEN / DRIED |
| 5. Variety | DIFFERENT THAN CALIBRATORS |

III. PRECAUTIONS:

1. Although the urine has been tested and found negative for HBsAg by RIA and HIV by EIA, the control material should be treated as any other potentially infectious agent.
2. For *in vitro* diagnostic use only.

IV. STORAGE AND STABILITY:

1. Store dried control material at 2-8°C (35-46°F). Stable to expiration date printed on the insert and label.
2. Store reconstituted control material at or below 2-8°C (35-46°F). See table for reconstituted stability.

V. PROCEDURE:

1. Remove cap from each vial to be used.
2. Reconstitute control material by adding exactly **5 mL** of **distilled water**, using a **5 mL** volumetric pipette or equivalent.
3. Replace cap and let sit 10-15 minutes.
4. Swirl gently 3-4 minutes to ensure a homogeneous mixture.
5. Swirl gently each time an aliquot is removed to ensure a homogeneous mixture.

VI. LIMITATIONS:

1. Results are dependent upon proper storage and adequate mixing.
2. Control material approximates a patient specimen.

VII. EXPECTED VALUES:

1. Laboratories should establish their own mean values; an individual laboratory's mean of several determinations may not duplicate the values listed below.

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C8582

Lot Number:

Product Description: Benzodiazepines Plus 100 ng/mL

2024-03-31

Expiration Date:

Analyte	Stability(Days)	Gravimetrically Weighed Value	Verified Value	Units
Alprazolam	25	100	95	ng/mL
αHydroxyalprazolam	25	100	89	ng/mL
Chlordiazepoxide	25	100	87	ng/mL
Norchlordiazepoxide	25	100	Present	ng/mL
Clobazam	25	100	100	ng/mL
Clonazepam	25	100	100	ng/mL
7-Aminoclonazepam	N/A*	100	81	ng/mL
Diazepam	25	100	110	ng/mL
Nordiazepam	25	100	95	ng/mL
Estazolam	25	100	100	ng/mL
Flunitrazepam	25	100	87	ng/mL
7-Aminoflunitrazepam	N/A*	100	98	ng/mL
Flurazepam	25	100	110	ng/mL
Desalkylflurazepam	25	100	95	ng/mL
Hydroxyethylflurazepam	25	100	100	ng/mL
Lorazepam	25	100	98	ng/mL
Midazolam	25	100	102	ng/mL
αHydroxymidazolam	25	100	91	ng/mL
Nitrazepam	25	100	100	ng/mL
7-Aminonitrazepam	N/A*	100	100	ng/mL
Oxazepam	25	100	97	ng/mL
Temazepam	25	100	93	ng/mL
Triazolam	25	100	94	ng/mL
αHydroxytriazolam	25	100	100	ng/mL
Zolpidem	25	100	100	ng/mL

Matrix: Urine

Product Form: Lyophilized

Storage: 2-8°C (35-46°F)

Vial Volume: 5 mL

*7-Amino metabolites will be weighed and spiked into this control material; however, no claim will be made for expected values or stability.

UTAK's express and implied warranties (including merchantability and fitness) are conditioned on the observance of UTAK's insert directions with respect to the use of UTAK's products.

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For technical support please call: UTAK Technical Service (888) 882-5522

IVD

UTAK Laboratories, Inc.
25020 Avenue Tibbitts
Valencia, CA 91355
T: (888) 882-5522
F: (661) 294-9272
E: welovecontrol@utak.com

EMERGO EUROPE
PRINSESSEGRACHT 20
2514 AP THE HAGUE
THE NETHERLANDS