ENGLISH:

I. INTENDED USE:

The Benzodiazepines Plus 100 ng/mL, urine toxicology control is for use as a control material to evaluate the performance of an analytical test method.

II. SUMMARY AND PRINCIPLES:

Several different techniques are used for evaluating or estimating the variance of results. The three subjects summarized below must be considered with any test method.

1. PREVENTIVE MEASURES:

These measures are usually contained in the design of the test method and include consideration for reagents, equipment, and operator errors. These measures are designed to minimize variance.

2. CONTROL MEASURES:

When a control sample is analyzed at the same time and in the same manner as a patient specimen, an estimate of variance is obtained for the test method. This estimate of variance can be compared to the acceptable limits of variance of the test method.

3. STATISTICAL ANALYSIS OF PATIENT RESULTS:

As an aid in evaluating overall test results, the past experience of expected results can be compared to the results of any given test run. For example, it would not be expected that all results of a given test run be in an elevated range.

Control materials are widely used as a means to aid in the evaluation of test results. The following subjects are to be considered in the use of any control material.

1. Multi-Level	NORMAL / ELEVATED
2. Matrix	HUMAN / ANIMAL / CHEMICAL
3. Availability	SUFFICIENT FOR STATISTICS
4. Form	LIQUID / FROZEN / DRIED
5. Variety	DIFFERENT THAN CALIBRATORS

III. PRECAUTIONS:

- Although the urine has been tested and found negative for HBsAg by RIA and HIV by EIA, the control material should be treated as any other potentially infectious agent.
- For *in vitro* diagnostic use only.

IV. STORAGE AND STABILITY:

- Store dried control material at 2-8°C (35-46°F). Stable to expiration date printed on the insert and label.
- Store reconstituted control material at or below 2-8°C (35-46°F). See table for reconstituted stability.

V. PROCEDURE:

- Remove cap from each vial to be used.
- Reconstitute control material by adding exactly 5 mL of distilled water, using a 5 mL volumetric pipette or equivalent.
- Replace cap and let sit 10-15 minutes.
- Swirl gently 3-4 minutes to ensure a homogeneous mixture.
- Swirl gently each time an aliquot is removed to ensure a homogeneous mixture.

VI. LIMITATIONS:

- Results are dependent upon proper storage and adequate mixing.
- Control material approximates a patient specimen.

VII. EXPECTED VALUES:

- Laboratories should establish their own mean values; an individual laboratory's mean of several determinations may not duplicate the values listed below.

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REF 12090 Product Number:		Product Description: Benzodiazepines Plus 100 ng/mL		
LOT C9839 Lot Number:		 2024-09-30 Expiration Date:		
Analyte	Stability(Days)	Gravimetrically Weighed Value	Verified Value	Units
Alprazolam	25	100	96	ng/mL
α-Hydroxyalprazolam	25	100	84.5	ng/mL
Chlordiazepoxide	25	100	84	ng/mL
Norchlordiazepoxide	25	100	Present	ng/mL
Clobazam	25	100	94.5	ng/mL
Clonazepam	25	100	87	ng/mL
7-Aminoclonazepam	N/A*	100	85	ng/mL
Diazepam	25	100	96	ng/mL
Nordiazepam	25	100	97	ng/mL
Estazolam	25	100	110	ng/mL
Flunitrazepam	25	100	93	ng/mL
7-Aminoflunitrazepam	N/A*	100	95.5	ng/mL
Flurazepam	25	100	95	ng/mL
Desalkylflurazepam	25	100	100	ng/mL
Hydroxyethylflurazepam	25	100	86	ng/mL
Lorazepam	25	100	93	ng/mL
Midazolam	25	100	89	ng/mL
α-Hydroxymidazolam	25	100	90.5	ng/mL
Nitrazepam	25	100	100	ng/mL
7-Aminonitrazepam	N/A*	100	86	ng/mL
Oxazepam	25	100	94	ng/mL
Temazepam	25	100	95	ng/mL
Triazolam	25	100	95	ng/mL
α-Hydroxytriazolam	25	100	94.5	ng/mL
Zolpidem	25	100	103	ng/mL
Matrix: Urine Product Form: Lyophilized Storage: 2-8°C (35-46°F) Vial Volume: 5 mL				

*7-Amino metabolites will be weighed and spiked into this control material; however, no claim will be made for expected values or stability.

UTAK's express and implied warranties (including merchantability and fitness) are conditioned on the observance of UTAK's insert directions with respect to the use of UTAK's products.

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For technical support please call: UTAK Technical Service (888) 882-5522

IVD

UTAK Laboratories, Inc.
25020 Avenue Tibbitts
Valencia, CA 91355
T: (888) 882-5522
F: (661) 294-9272
E: welovecontrol@utak.com

EMERGO EUROPE
PRINSESSEGRACHT 20
2514 AP THE HAGUE
THE NETHERLANDS