REF 12090

Product Description: Benzodiazepines Plus 100
ng/mL

Product Number:

LOT D7594

Lot Number:

 2028-03-31

Expiration Date:

Analyte	Reconstituted Stability (Days)	Gravimetrically Weighed Value	Verified Value	Units
Alprazolam	25	100	93.5	ng/mL
α-Hydroxyalprazolam	25	100	96	ng/mL
Chlordiazepoxide	25	100	88.5	ng/mL
Norchlordiazepoxide	25	100	Present	ng/mL
Clobazam	25	100	100	ng/mL
Clonazepam	25	100	91	ng/mL
7-Aminoclonazepam	N/A ¹	100	63.3	ng/mL
Diazepam	25	100	94.5	ng/mL
Nordiazepam	25	100	97	ng/mL
Estazolam	25	100	100	ng/mL
Flunitrazepam	25	100	96	ng/mL
7-Aminoflunitrazepam	N/A ¹	100	83.5	ng/mL
Flurazepam	25	100	103	ng/mL
Desalkylflurazepam	25	100	100	ng/mL
Hydroxyethylflurazepam	25	100	100.5	ng/mL
Lorazepam	25	100	92.5	ng/mL
Midazolam	25	100	98	ng/mL
α-Hydroxymidazolam	25	100	98	ng/mL
Nitrazepam	25	100	100	ng/mL
7-Aminonitrazepam	N/A ¹	100	61	ng/mL
Oxazepam	25	100	101	ng/mL
Temazepam	25	100	100.5	ng/mL
Triazolam	25	100	96	ng/mL
α-Hydroxytriazolam	25	100	97.5	ng/mL
Zolpidem	25	100	103	ng/mL

Matrix: Urine

Product Form: Lyophilized

Storage: 2-8°C (35-46°F)

Vial Volume: 5 mL

¹ Analyte has been weighed and spiked into the control; however, no claim will be made for expected values or stability.

ENGLISH:

I. INTENDED USE:

The 12090, Benzodiazepines Plus 100 ng/mL, urine control is for use as a control is for use as a control material to evaluate the performance of an analytical test method.

II. SUMMARY AND PRINCIPLES:

Several different techniques are used for evaluating or estimating the variance of results. The three subjects summarized below must be considered with any test method.

1. PREVENTIVE MEASURES:

These measures are usually contained in the design of the test method and include consideration for reagents, equipment, and operator errors. These measures are designed to minimize variance.

2. CONTROL MEASURES:

When a control sample is analyzed at the same time and in the same manner as a patient specimen, an estimate of variance is obtained for the test method. This estimate of variance can be compared to the acceptable limits of variance of the test method.

3. STATISTICAL ANALYSIS OF PATIENT RESULTS:

As an aid in evaluating overall test results, the past experience of expected results can be compared to the results of any given test run. For example, it would not be expected that all results of a given test run be in an elevated range.

Control materials are widely used as a means to aid in the evaluation of test results. The following subjects are to be considered in the use of any control material.

1. Multi-Level	NORMAL / ELEVATED
2. Matrix	HUMAN / ANIMAL / CHEMICAL
3. Availability	SUFFICIENT FOR STATISTICS
4. Form	LIQUID / FROZEN / DRIED
5. Variety	DIFFERENT THAN CALIBRATORS

III. PRECAUTIONS:

1. Although the urine has been tested and found negative for HBsAg by RIA and HIV by EIA, the control material should be treated as any other potentially infectious agent.
2. For *in vitro* diagnostic use only.

IV. STORAGE AND STABILITY:

1. Store dried control material at 2-8°C (35-46°F). Stable to expiration date printed on the insert and label.
2. Store reconstituted control material at 2-8°C (35-46°F). See table for reconstituted stability.

V. PROCEDURE:

1. Remove cap from each vial to be used.
2. Reconstitute control material by adding exactly **5 mL** of **distilled water**, using a **5 mL** volumetric pipette or equivalent.
3. Replace cap and let sit 10-15 minutes.
4. Swirl gently 3-4 minutes to ensure a homogeneous mixture.
5. Swirl gently each time an aliquot is removed to ensure a homogeneous mixture.

VI. LIMITATIONS:

1. Results are dependent upon proper storage and adequate mixing.
2. Control material approximates a patient specimen.

VII. EXPECTED VALUES:

1. Laboratories should establish their own mean values; an individual laboratory's mean of several determinations may not duplicate the values listed below.

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