### ENGLISH:

#### I. INTENDED USE:

The Benzodiazepines Plus 100 ng/mL, whole blood toxicology control is for use as a control material to evaluate the performance of an analytical test method.

#### II. SUMMARY AND PRINCIPLES:

Several different techniques are used for evaluating or estimating the variance of results. The three subjects summarized below must be considered with any test method.

##### 1. PREVENTIVE MEASURES:

These measures are usually contained in the design of the test method and include consideration for reagents, equipment, and operator errors. These measures are designed to minimize variance.

##### 2. CONTROL MEASURES:

When a control sample is analyzed at the same time and in the same manner as a patient specimen, an estimate of variance is obtained for the test method. This estimate of variance can be compared to the acceptable limits of variance of the test method.

##### 3. STATISTICAL ANALYSIS OF PATIENT RESULTS:

As an aid in evaluating overall test results, the past experience of expected results can be compared to the results of any given test run. For example, it would not be expected that all results of a given test run be in an elevated range.

Control materials are widely used as a means to aid in the evaluation of test results. The following subjects are to be considered in the use of any control material.

1. Multi-Level	NORMAL / ELEVATED
2. Matrix	HUMAN / ANIMAL / CHEMICAL
3. Availability	SUFFICIENT FOR STATISTICS
4. Form	LIQUID / FROZEN / DRIED
5. Variety	DIFFERENT THAN CALIBRATORS

#### III. PRECAUTIONS:

- Although the whole blood has been tested and found negative for HBsAg by RIA and HIV by EIA, the control material should be treated as any other potentially infectious agent.
- For *in vitro* diagnostic use only.

#### IV. STORAGE AND STABILITY:

- Store dried control material at 2-8°C (35-46°F). Stable to expiration date printed on the insert and label.
- Store reconstituted control material at or below 2-8°C (35-46°F). See table for reconstituted stability.

#### V. PROCEDURE:

- Remove cap from each vial to be used.
- Reconstitute control material by adding exactly 5 mL of distilled water, using a 5 mL volumetric pipette or equivalent.
- Replace cap and let sit 10-15 minutes.
- Swirl gently 3-4 minutes to ensure a homogeneous mixture.
- Swirl gently each time an aliquot is removed to ensure a homogeneous mixture.

#### VI. LIMITATIONS:

- Results are dependent upon proper storage and adequate mixing.
- Control material approximates a patient specimen.

#### VII. EXPECTED VALUES:

- Laboratories should establish their own mean values; an individual laboratory's mean of several determinations may not duplicate the values listed below.

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<b>REF</b>	12100	<b>Product Description:</b> Benzodiazepines Plus 100 ng/mL
Product Number:		
<b>LOT</b>	D1949	2025-03-31
Lot Number:		Expiration Date:

Analyte	Stability(Days)	Gravimetrically Weighed Value	Verified Value	Units
Alprazolam	25	100	76	ng/mL
α-Hydroxyalprazolam	25	100	85	ng/mL
Chlordiazepoxide	25	100	76	ng/mL
Clobazam	25	100	66	ng/mL
Clonazepam	25	100	71	ng/mL
7-Aminoclonazepam	N/A*	100	37	ng/mL
Diazepam	25	100	89	ng/mL
Nordiazepam	25	100	71	ng/mL
Estazolam	25	100	80	ng/mL
Flunitrazepam	25	100	35	ng/mL
7-Aminoflunitrazepam	N/A*	100	66	ng/mL
Flurazepam	25	100	67	ng/mL
Desalkylflurazepam	25	100	79	ng/mL
2-Hydroxyethylflurazepam	25	100	67	ng/mL
Lorazepam	25	100	79	ng/mL
Midazolam	25	100	85	ng/mL
Nitrazepam	25	100	80	ng/mL
7-Aminonitrazepam	N/A*	100	36	ng/mL
Oxazepam	25	100	74	ng/mL
Temazepam	25	100	68	ng/mL
Triazolam	25	100	76	ng/mL
α-Hydroxytriazolam	25	100	84	ng/mL
Zolpidem	25	100	68	ng/mL

**Matrix:** Whole Blood  
**Product Form:** Lyophilized  
**Storage:** 2-8°C (35-46°F)  
**Vial Volume:** 5 mL

\*7-Amino metabolites have been weighed and spiked into the control; however, no claim will be made for expected values or stability.

UTAK's express and implied warranties (including merchantability and fitness) are conditioned on the observance of UTAK's insert directions with respect to the use of UTAK's products.

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For technical support please call: UTAK Technical Service (888) 882-5522

IVD

UTAK Laboratories, Inc.  
25020 Avenue Tibbitts  
Valencia, CA 91355  
T: (888) 882-5522  
F: (661) 294-9272  
E: [welovecontrol@utak.com](mailto:welovecontrol@utak.com)

EMERGO EUROPE  
PRINSESSEGRACHT 20  
2514 AP THE HAGUE  
THE NETHERLANDS