REF 12110

Product Description: HR Metals

Product Number:

LOT D5343

Lot Number:

2027-01-31

Expiration Date:

Analyte	Reconstituted Stability (Days)	Method	Verified Value	Expected Range	Units
Aluminum	25	ICP/MS	46	39.1 - 52.9	µg/L
Arsenic	25	ICP/MS	113	96.05 - 129.95	µg/L
Cadmium	25	ICP/MS	4.9	4.17 - 5.64	µg/L
Calcium	25	ICP/OES	433	368.05 - 497.95	mg/L
Chromium	25	ICP/MS	2.5	2.13 - 2.88	µg/L
Cobalt	25	ICP/MS	6.2	5.27 - 7.13	µg/L
Copper	25	ICP/MS	77	65.45 - 88.55	µg/L
Creatinine	25	JAFFE	1.663	1.41 - 1.91	g/L
Fluoride	25	ISE	790	671.5 - 908.5	µg/L
Iron	25	ICP/OES	660	561 - 759	µg/L
Lead	25	ICP/MS	184	156.4 - 211.6	µg/L
Manganese	25	ICP/MS	1.1	0.94 - 1.27	µg/L
Mercury	25	ICP/MS	40	34 - 46	µg/L
Molybdenum	25	ICP/MS	67	56.95 - 77.05	µg/L
Nickel	25	ICP/MS	30	25.5 - 34.5	µg/L
Selenium	25	ICP/MS	46	39.1 - 52.9	µg/L
Vanadium	25	ICP/MS	12	10.2 - 13.8	µg/L
Zinc	25	ICP/MS	970	824.5 - 1115.5	µg/L

Matrix: Urine

Product Form: Lyophilized

Storage: 2-8°C (35-46°F)

Vial Volume: 5 mL

ENGLISH:

I. INTENDED USE:

The 12110, HR Metals, urine control is for use as a control material to evaluate the performance of an analytical test method.

II. SUMMARY AND PRINCIPLES:

Several different techniques are used for evaluating or estimating the variance of results. The three subjects summarized below must be considered with any test method.

1. PREVENTIVE MEASURES:

These measures are usually contained in the design of the test method and include consideration for reagents, equipment, and operator errors. These measures are designed to minimize variance.

2. CONTROL MEASURES:

When a control sample is analyzed at the same time and in the same manner as a patient specimen, an estimate of variance is obtained for the test method. This estimate of variance can be compared to the acceptable limits of variance of the test method.

3. STATISTICAL ANALYSIS OF PATIENT RESULTS:

As an aid in evaluating overall test results, the past experience of expected results can be compared to the results of any given test run. For example, it would not be expected that all results of a given test run be in an elevated range.

Control materials are widely used as a means to aid in the evaluation of test results. The following subjects are to be considered in the use of any control material.

- | | |
|-----------------|----------------------------|
| 1. Multi-Level | NORMAL / ELEVATED |
| 2. Matrix | HUMAN / ANIMAL / CHEMICAL |
| 3. Availability | SUFFICIENT FOR STATISTICS |
| 4. Form | LIQUID / FROZEN / DRIED |
| 5. Variety | DIFFERENT THAN CALIBRATORS |

III. PRECAUTIONS:

- Although the urine has been tested and found negative for HBsAg by RIA and HIV by EIA, the control material should be treated as any other potentially infectious agent.
- For *in vitro* diagnostic use only.

IV. STORAGE AND STABILITY:

- Store dried control material at 2-8°C (35-46°F). Stable to expiration date printed on the insert and label.
- Store reconstituted control material at 2-8°C (35-46°F). See table for reconstituted stability.

V. PROCEDURE:

- Remove cap from each vial to be used.
- Reconstitute control material by adding exactly **5 mL** of **distilled water**, using a **5 mL** volumetric pipette or equivalent.
- Replace cap and let sit 10-15 minutes.
- Swirl gently 3-4 minutes to ensure a homogeneous mixture.
- Swirl gently each time an aliquot is removed to ensure a homogeneous mixture.

VI. LIMITATIONS:

- Results are dependent upon proper storage and adequate mixing.
- Control material approximates a patient specimen.

VII. EXPECTED VALUES:

- Laboratories should establish their own mean values; an individual laboratory's mean of several determinations may not duplicate the values listed below.

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