ENGLISH:**I. INTENDED USE:**

The Validity Control 5, adulterant toxicology control is for use as a control material to evaluate the performance of an analytical test method.

II. SUMMARY AND PRINCIPLES:

Several different techniques are used for evaluating or estimating the variance of results. The three subjects summarized below must be considered with any test method.

1. PREVENTIVE MEASURES:

These measures are usually contained in the design of the test method and include consideration for reagents, equipment, and operator errors. These measures are designed to minimize variance.

2. CONTROL MEASURES:

When a control sample is analyzed at the same time and in the same manner as a patient specimen, an estimate of variance is obtained for the test method. This estimate of variance can be compared to the acceptable limits of variance of the test method.

3. STATISTICAL ANALYSIS OF PATIENT RESULTS:

As an aid in evaluating overall test results, the past experience of expected results can be compared to the results of any given test run. For example, it would not be expected that all results of a given test run be in an elevated range.

Control materials are widely used as a means to aid in the evaluation of test results. The following subjects are to be considered in the use of any control material.

1. Multi-Level	NORMAL / ELEVATED
2. Matrix	HUMAN / ANIMAL / CHEMICAL
3. Availability	SUFFICIENT FOR STATISTICS
4. Form	LIQUID / FROZEN / DRIED
5. Variety	DIFFERENT THAN CALIBRATORS

III. PRECAUTIONS:**1. For *in vitro* diagnostic use only.****IV. STORAGE AND STABILITY:****1. Store liquid control material at 2-8°C (35-46°F). Stable to expiration date printed on the label.****V. PROCEDURE:**

1. Swirl gently 3-4 minutes to ensure a homogeneous mixture.
2. Swirl gently each time an aliquot is removed to ensure a homogeneous mixture.
3. Assay control material in same manner as patient specimens, following the exact same instructions from the entire test method.
4. Record the results obtained on a control chart that describes statistical limits for the test method and the particular lot of control material.

VI. LIMITATIONS:

1. Results are dependent upon proper storage and adequate mixing.
2. Control material approximates a patient specimen.

VII. EXPECTED VALUES:

1. Laboratories should establish their own mean values; an individual laboratory's mean of several determinations may not duplicate the values listed below.

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REF

17014

Product Description: Validity Control 5

Product Number:

LOT

D1144

Lot Number:

2024-05-31

Expiration Date:

Analyte	Stability(Days)	Expected Range	Units
Chromium VI	to expiration date	>65	µg/mL
Creatinine	to expiration date	n/a	mg/dL
Nitrites	to expiration date	n/a	µg/mL
Oxidants	to expiration date	>65	µg/mL
pH	to expiration date	11.2 - 12.0	n/a
Specific Gravity	to expiration date	n/a	n/a
Matrix: Adulterant Product Form: Liquid Storage: 2-8°C (35-46°F) Vial Volume: 25 mL			

UTAK's express and implied warranties (including merchantability and fitness) are conditioned on the observance of UTAK's insert directions with respect to the use of UTAK's products.

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For technical support please call: UTAK Technical Service (888) 882-5522

IVD

UTAK Laboratories, Inc.
25020 Avenue Tibbitts
Valencia, CA 91355
T: (888) 882-5522
F: (661) 294-9272
E: welovecontrol@utak.com

EMERGO EUROPE
PRINSESSEGRACHT 20
2514 AP THE HAGUE
THE NETHERLANDS