ENGLISH:

I. INTENDED USE:

The Amphetamines Panel Level 2, urine toxicology control is for use as a control material to evaluate the performance of an analytical test method.

II. SUMMARY AND PRINCIPLES:

Several different techniques are used for evaluating or estimating the variance of results. The three subjects summarized below must be considered with any test method.

1. PREVENTIVE MEASURES:

These measures are usually contained in the design of the test method and include consideration for reagents, equipment, and operator errors. These measures are designed to minimize variance.

2. CONTROL MEASURES:

When a control sample is analyzed at the same time and in the same manner as a patient specimen, an estimate of variance is obtained for the test method. This estimate of variance can be compared to the acceptable limits of variance of the test method.

3. STATISTICAL ANALYSIS OF PATIENT RESULTS:

As an aid in evaluating overall test results, the past experience of expected results can be compared to the results of any given test run. For example, it would not be expected that all results of a given test run be in an elevated range.

Control materials are widely used as a means to aid in the evaluation of test results. The following subjects are to be considered in the use of any control material.

1. Multi-Level	NORMAL / ELEVATED
2. Matrix	HUMAN / ANIMAL / CHEMICAL
3. Availability	SUFFICIENT FOR STATISTICS
4. Form	LIQUID / FROZEN / DRIED
5. Variety	DIFFERENT THAN CALIBRATORS

III. PRECAUTIONS:

1. Although the urine has been tested and found negative for HBsAg by RIA and HIV by EIA, the control material should be treated as any other potentially infectious agent.
2. For *in vitro* diagnostic use only.

IV. STORAGE AND STABILITY:

1. Store dried control material at 2-8°C (35-46°F). Stable to expiration date printed on the insert and label.
2. Store reconstituted control material at or below 2-8°C (35-46°F). See table for reconstituted stability.

V. PROCEDURE:

1. Remove cap from each vial to be used.
2. Reconstitute control material by adding exactly **5 mL** of **distilled water**, using a **5 mL** volumetric pipette or equivalent.
3. Replace cap and let sit 10-15 minutes.
4. Swirl gently 3-4 minutes to ensure a homogeneous mixture.
5. Swirl gently each time an aliquot is removed to ensure a homogeneous mixture.

VI. LIMITATIONS:

1. Results are dependent upon proper storage and adequate mixing.
2. Control material approximates a patient specimen.

VII. EXPECTED VALUES:

1. Laboratories should establish their own mean values; an individual laboratory's mean of several determinations may not duplicate the values listed below.

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REF	41401	Product Description: Amphetamines Panel Level 2
Product Number:		
LOT	D3392	 2025-12-31 Expiration Date:
Lot Number:		

Analyte	Stability(Days)	Gravimetrically Weighed Value	Verified Value	Units
d-Amphetamine	25	700	620	ng/mL
(±)-3,4-Methylenedioxyamphetamine (MDA)	25	700	670	ng/mL
d-Methamphetamine	25	700	690	ng/mL
(±)-3,4-Methylenedioxymethamphetamine (MDMA)	25	700	680	ng/mL
(±)-3,4-Methylenedioxyethylamphetamine (MDEA)	25	700	660	ng/mL
(-)-Ephedrine	25	700	710	ng/mL
Methylephedrine	25	700	Present	ng/mL
Phendimetrazine	25	700	Present	ng/mL
Phenmetrazine	25	700	Present	ng/mL
Phentermine	25	700	660	ng/mL
Phenylpropanolamine	25	700	590	ng/mL
(+)-Pseudoephedrine	25	700	640	ng/mL
Selegiline	25	700	Present	ng/mL

Matrix: Urine
Product Form: Lyophilized
Storage: 2-8°C (35-46°F)
Vial Volume: 5 mL

UTAK's express and implied warranties (including merchantability and fitness) are conditioned on the observance of UTAK's insert directions with respect to the use of UTAK's products.

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