### ENGLISH:

#### I. INTENDED USE:

The Drugs Of Abuse 50% Cutoff, oral fluid toxicology control is for use as a control material to evaluate the performance of an analytical test method.

#### II. SUMMARY AND PRINCIPLES:

Several different techniques are used for evaluating or estimating the variance of results. The three subjects summarized below must be considered with any test method.

##### 1. PREVENTIVE MEASURES:

These measures are usually contained in the design of the test method and include consideration for reagents, equipment, and operator errors. These measures are designed to minimize variance.

##### 2. CONTROL MEASURES:

When a control sample is analyzed at the same time and in the same manner as a patient specimen, an estimate of variance is obtained for the test method. This estimate of variance can be compared to the acceptable limits of variance of the test method.

##### 3. STATISTICAL ANALYSIS OF PATIENT RESULTS:

As an aid in evaluating overall test results, the past experience of expected results can be compared to the results of any given test run. For example, it would not be expected that all results of a given test run be in an elevated range.

Control materials are widely used as a means to aid in the evaluation of test results. The following subjects are to be considered in the use of any control material.

|                 |                            |
|-----------------|----------------------------|
| 1. Multi-Level  | NORMAL / ELEVATED          |
| 2. Matrix       | HUMAN / ANIMAL / CHEMICAL  |
| 3. Availability | SUFFICIENT FOR STATISTICS  |
| 4. Form         | LIQUID / FROZEN / DRIED    |
| 5. Variety      | DIFFERENT THAN CALIBRATORS |

#### III. PRECAUTIONS:

1. For *in vitro* diagnostic use only.

#### IV. STORAGE AND STABILITY:

1. Store fresh frozen control material at or below -10°C (14°F). Stable to expiration date printed on the insert and label.
2. Store thawed control material at or below 2-8°C (35-46°F). See table for thawed stability.

#### V. PROCEDURE:

1. Allow fresh frozen control material to thaw at room temperature with cap on.
2. Swirl gently 3-4 minutes to ensure a homogeneous mixture.
3. Swirl gently each time an aliquot is removed to ensure a homogeneous mixture.
4. Assay control material in same manner as patient specimens, following the exact same instructions from the entire test method.
5. Record the results obtained on a control chart that describes statistical limits for the test method and the particular lot of control material.

#### VI. LIMITATIONS:

1. Results are dependent upon proper storage and adequate mixing.
2. Control material approximates a patient specimen.

#### VII. EXPECTED VALUES:

1. Laboratories should establish their own mean values; an individual laboratory's mean of several determinations may not duplicate the values listed below.

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|                        |       |   |
|------------------------|-------|---|
| <b>REF</b>             | 43050 | <b>Product Description:</b> Drugs Of Abuse 50% Cutoff |
| <b>Product Number:</b> |       |   |
| <b>LOT</b>             | D1179 | 2023-05-31  |
| <b>Lot Number:</b>     |       | <b>Expiration Date:</b>                               |

| Analyte             | Thawed Stability<br>(Days) | Gravimetrically<br>Weighed Value | Verified Value | Units |
|---------------------|----------------------------|----------------------------------|----------------|-------|
| d-Amphetamine       | 25                         | 25                               | 26.9           | ng/mL |
| d-Methamphetamine   | 25                         | 25                               | 27             | ng/mL |
| Benzoylcegonine     | 25                         | 10                               | 14.85          | ng/mL |
| Methadone           | 25                         | 2.5                              | 3              | ng/mL |
| Morphine            | 25                         | 20                               | 21.6           | ng/mL |
| Oxazepam            | 25                         | 25                               | 28.7           | ng/mL |
| Phencyclidine (PCP) | 25                         | 5                                | 6              | ng/mL |
| Propoxyphene        | 25                         | 20                               | 20.4           | ng/mL |
| Secobarbital        | 25                         | 15                               | 15.7           | ng/mL |
| Delta-9-THC         | 14                         | 2                                | 2.75           | ng/mL |

**Matrix:** Oral Fluid  
**Product Form:** Frozen  
**Storage:** -10°C (14°F)  
**Vial Volume:** 3 mL

UTAK's express and implied warranties (including merchantability and fitness) are conditioned on the observance of UTAK's insert directions with respect to the use of UTAK's products.

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