REF 50700

Product Description: DAU High Cutoff 1

Product Number:

LOT D5370

Lot Number:

 2027-05-31

Expiration Date:

Analyte	Reconstituted Stability (Days)	Gravimetrically Weighed Value	Verified Value	Units
α-Hydroxyalprazolam	25	225	237	ng/mL
Amobarbital	25	225	244	ng/mL
d-Amphetamine	25	375	406	ng/mL
(±)-3,4-Methylenedioxyamphetamine (MDA)	25	375	347	ng/mL
d-Methamphetamine	25	375	387	ng/mL
(±)-3,4-Methylenedioxymethamphetamine (MDMA)	25	375	378	ng/mL
(±)-3,4-Methylenedioxyethylamphetamine (MDEA)	25	375	354	ng/mL
Benzoyllecgonine	25	112.5	100	ng/mL
Butalbital	25	225	237	ng/mL
Codeine	25	1500	1343	ng/mL
Hydrocodone	25	300	242	ng/mL
Hydromorphone	25	300	280	ng/mL
Nordiazepam	25	225	256	ng/mL
Lorazepam	25	225	246	ng/mL
Methadone	25	225	239	ng/mL
Methadone Primary Metabolite (EDDP)	25	225	230	ng/mL
Methaqualone	25	225	219	ng/mL
Morphine-3-β-D-Glucuronide	25	2430 (Equivalent to 1500 free Morphine)	1346.22 (Equivalent to 831 free Morphine)	ng/mL
6-Monoacetylmorphine (6-MAM)	25	7.5	9	ng/mL
Nortriptyline	25	225	192	ng/mL
Oxazepam	25	225	236	ng/mL
Oxycodone	25	225	206	ng/mL
Oxymorphone	25	225	206	ng/mL
Pentobarbital	25	225	224	ng/mL
Phencyclidine (PCP)	25	18.75	20	ng/mL
Phenobarbital	25	225	278	ng/mL
Norpropoxyphene	25	225	237	ng/mL
Secobarbital	25	225	215	ng/mL
Temazepam	25	225	238	ng/mL
(l)-9-Carboxy-11-nor-Delta-9-THC (THC-COOH)	25	11.25	9	ng/mL
Matrix: Urine Product Form: Lyophilized Storage: 2-8°C (35-46°F) Vial Volume: 10 mL				

ENGLISH:

I. INTENDED USE:

The 50700, DAU High Cutoff 1, urine control is for use as a control material to evaluate the performance of an analytical test method.

II. SUMMARY AND PRINCIPLES:

Several different techniques are used for evaluating or estimating the variance of results. The three subjects summarized below must be considered with any test method.

1. PREVENTIVE MEASURES:

These measures are usually contained in the design of the test method and include consideration for reagents, equipment, and operator errors. These measures are designed to minimize variance.

2. CONTROL MEASURES:

When a control sample is analyzed at the same time and in the same manner as a patient specimen, an estimate of variance is obtained for the test method. This estimate of variance can be compared to the acceptable limits of variance of the test method.

3. STATISTICAL ANALYSIS OF PATIENT RESULTS:

As an aid in evaluating overall test results, the past experience of expected results can be compared to the results of any given test run. For example, it would not be expected that all results of a given test run be in an elevated range.

Control materials are widely used as a means to aid in the evaluation of test results. The following subjects are to be considered in the use of any control material.

1. Multi-Level	NORMAL / ELEVATED
2. Matrix	HUMAN / ANIMAL / CHEMICAL
3. Availability	SUFFICIENT FOR STATISTICS
4. Form	LIQUID / FROZEN / DRIED
5. Variety	DIFFERENT THAN CALIBRATORS

III. PRECAUTIONS:

1. Although the urine has been tested and found negative for HBsAg by RIA and HIV by EIA, the control material should be treated as any other potentially infectious agent.
2. For *in vitro* diagnostic use only.

IV. STORAGE AND STABILITY:

1. Store dried control material at 2-8°C (35-46°F). Stable to expiration date printed on the insert and label.
2. Store reconstituted control material at 2-8°C (35-46°F). See table for reconstituted stability.

V. PROCEDURE:

1. Remove cap from each vial to be used.
2. Reconstitute control material by adding exactly **10 mL of distilled water**, using a **10 mL** volumetric pipette or equivalent.
3. Replace cap and let sit 10-15 minutes.
4. Swirl gently 3-4 minutes to ensure a homogeneous mixture.
5. Swirl gently each time an aliquot is removed to ensure a homogeneous mixture.

VI. LIMITATIONS:

1. Results are dependent upon proper storage and adequate mixing.
2. Control material approximates a patient specimen.

VII. EXPECTED VALUES:

1. Laboratories should establish their own mean values; an individual laboratory's mean of several determinations may not duplicate the values listed below.

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