<b>REF</b>	50701	Product Description: DAU High Cutoff 2
Product Number:		
<b>LOT</b>	D8286	 2028-09-30 Expiration Date:
Lot Number:		

Analyte	Reconstituted Stability (Days)	Gravimetrically Weighed Value	Verified Value	Units
α-Hydroxyalprazolam	25	375	446	ng/mL
Amobarbital	25	375	432	ng/mL
d-Amphetamine	25	625	625	ng/mL
(±)-3,4-Methylenedioxyamphetamine (MDA)	25	625	597	ng/mL
d-Methamphetamine	25	625	651	ng/mL
(±)-3,4-Methylenedioxymethamphetamine (MDMA)	25	625	599	ng/mL
(±)-3,4-Methylenedioxyethylamphetamine (MDEA)	25	625	621	ng/mL
Benzoylcegonine	25	187.5	192	ng/mL
Butalbital	25	375	399	ng/mL
Codeine	25	2500	2567	ng/mL
Hydrocodone	25	500	489	ng/mL
Hydromorphone	25	500	511	ng/mL
Nordiazepam	25	375	422	ng/mL
Lorazepam	25	375	406	ng/mL
Methadone	25	375	355	ng/mL
Methadone Primary Metabolite (EDDP)	25	375	391	ng/mL
Methaqualone	25	375	401	ng/mL
Morphine-3-β-D-Glucuronide	25	4050 (Equivalent to 2500 free Morphine)	3284 (Equivalent to 2027 free Morphine)	ng/mL
6-Monoacetylmorphine (6-MAM)	25	12.5	13	ng/mL
Nortriptyline	25	375	293	ng/mL
Oxazepam	25	375	442	ng/mL
Oxycodone	25	375	396	ng/mL
Oxymorphone	25	375	444	ng/mL
Pentobarbital	25	375	412	ng/mL
Phencyclidine (PCP)	25	31.25	33	ng/mL
Phenobarbital	25	375	385	ng/mL
Norpropoxyphene	25	375	380	ng/mL
Secobarbital	25	375	401	ng/mL
Temazepam	25	375	440	ng/mL
(l)-9-Carboxy-11-nor-Delta-9-THC (THC-COOH)	25	18.75	18	ng/mL

**Matrix:** Urine  
**Product Form:** Lyophilized  
**Storage:** 2-8°C (35-46°F)  
**Vial Volume:** 10 mL

### ENGLISH:

#### I. INTENDED USE:

The 50701, DAU High Cutoff 2, urine control is for use as a control material to evaluate the performance of an analytical test method.

#### II. SUMMARY AND PRINCIPLES:

Several different techniques are used for evaluating or estimating the variance of results. The three subjects summarized below must be considered with any test method.

##### 1. PREVENTIVE MEASURES:

These measures are usually contained in the design of the test method and include consideration for reagents, equipment, and operator errors. These measures are designed to minimize variance.

##### 2. CONTROL MEASURES:

When a control sample is analyzed at the same time and in the same manner as a patient specimen, an estimate of variance is obtained for the test method. This estimate of variance can be compared to the acceptable limits of variance of the test method.

##### 3. STATISTICAL ANALYSIS OF PATIENT RESULTS:

As an aid in evaluating overall test results, the past experience of expected results can be compared to the results of any given test run. For example, it would not be expected that all results of a given test run be in an elevated range.

Control materials are widely used as a means to aid in the evaluation of test results. The following subjects are to be considered in the use of any control material.

1. Multi-Level	NORMAL / ELEVATED
2. Matrix	HUMAN / ANIMAL / CHEMICAL
3. Availability	SUFFICIENT FOR STATISTICS
4. Form	LIQUID / FROZEN / DRIED
5. Variety	DIFFERENT THAN CALIBRATORS

#### III. PRECAUTIONS:

- Although the urine has been tested and found negative for HBsAg by RIA and HIV by EIA, the control material should be treated as any other potentially infectious agent.
- For *in vitro* diagnostic use only.

#### IV. STORAGE AND STABILITY:

- Store dried control material at 2-8°C (35-46°F). Stable to expiration date printed on the insert and label.
- Store reconstituted control material at 2-8°C (35-46°F). See table for reconstituted stability.

#### V. PROCEDURE:

- Remove cap from each vial to be used.
- Reconstitute control material by adding exactly **10 mL of distilled water**, using a **10 mL** volumetric pipette or equivalent.
- Replace cap and let sit 10-15 minutes.
- Swirl gently 3-4 minutes to ensure a homogeneous mixture.
- Swirl gently each time an aliquot is removed to ensure a homogeneous mixture.

#### VI. LIMITATIONS:

- Results are dependent upon proper storage and adequate mixing.
- Control material approximates a patient specimen.

#### VII. EXPECTED VALUES:

- Laboratories should establish their own mean values; an individual laboratory's mean of several determinations may not duplicate the values listed below.

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**UTAK Laboratories, Inc.**