REF	50702	Product Description: DAU Low Cutoff 1
Product Number:		
LOT	D8597	2028-11-30 Expiration Date:
Lot Number:		

Analyte	Reconstituted Stability (Days)	Gravimetrically Weighed Value	Verified Value	Units
α-Hydroxyalprazolam	25	150	173	ng/mL
Amobarbital	25	150	130	ng/mL
d-Amphetamine	25	187.5	249	ng/mL
(±)-3,4-Methylenedioxyamphetamine (MDA)	25	187.5	209	ng/mL
d-Methamphetamine	25	187.5	222	ng/mL
(±)-3,4-Methylenedioxymethamphetamine (MDMA)	25	187.5	199	ng/mL
(±)-3,4-Methylenedioxyethylamphetamine (MDEA)	25	187.5	210	ng/mL
Benzoylcegonine	25	75	75	ng/mL
Butalbital	25	150	143	ng/mL
Codeine	25	225	223	ng/mL
Hydrocodone	25	75	74	ng/mL
Hydromorphone	25	75	71	ng/mL
Nordiazepam	25	150	167	ng/mL
Lorazepam	25	150	171	ng/mL
Methadone	25	150	156	ng/mL
Methadone Primary Metabolite (EDDP)	25	150	175	ng/mL
Methaqualone	25	150	170	ng/mL
Morphine-3-β-D-Glucuronide	25	365 (Equivalent to 225 free Morphine)	302.94 (Equivalent to 187 free Morphine)	ng/mL
6-Monoacetylmorphine (6-MAM)	25	7.5	8	ng/mL
Nortriptyline	25	150	146	ng/mL
Oxazepam	25	150	173	ng/mL
Oxycodone	25	37.5	36	ng/mL
Oxymorphone	25	37.5	37	ng/mL
Pentobarbital	25	150	139	ng/mL
Phencyclidine (PCP)	25	18.75	20	ng/mL
Phenobarbital	25	150	154	ng/mL
Norpropoxyphene	25	150	149	ng/mL
Secobarbital	25	150	147	ng/mL
Temazepam	25	150	159	ng/mL
(l)-9-Carboxy-11-nor-Delta-9-THC (THC-COOH)	25	11.25	10	ng/mL

Matrix: Urine
Product Form: Lyophilized
Storage: 2-8°C (35-46°F)
Vial Volume: 10 mL

ENGLISH:

I. INTENDED USE:

The 50702, DAU Low Cutoff 1, urine control is for use as a control material to evaluate the performance of an analytical test method.

II. SUMMARY AND PRINCIPLES:

Several different techniques are used for evaluating or estimating the variance of results. The three subjects summarized below must be considered with any test method.

- 1. PREVENTIVE MEASURES:**
These measures are usually contained in the design of the test method and include consideration for reagents, equipment, and operator errors. These measures are designed to minimize variance.
- 2. CONTROL MEASURES:**
When a control sample is analyzed at the same time and in the same manner as a patient specimen, an estimate of variance is obtained for the test method. This estimate of variance can be compared to the acceptable limits of variance of the test method.
- 3. STATISTICAL ANALYSIS OF PATIENT RESULTS:**
As an aid in evaluating overall test results, the past experience of expected results can be compared to the results of any given test run. For example, it would not be expected that all results of a given test run be in an elevated range.

Control materials are widely used as a means to aid in the evaluation of test results. The following subjects are to be considered in the use of any control material.

- | | |
|-----------------|----------------------------|
| 1. Multi-Level | NORMAL / ELEVATED |
| 2. Matrix | HUMAN / ANIMAL / CHEMICAL |
| 3. Availability | SUFFICIENT FOR STATISTICS |
| 4. Form | LIQUID / FROZEN / DRIED |
| 5. Variety | DIFFERENT THAN CALIBRATORS |

III. PRECAUTIONS:

1. Although the urine has been tested and found negative for HBsAg by RIA and HIV by EIA, the control material should be treated as any other potentially infectious agent.
2. For *in vitro* diagnostic use only.

IV. STORAGE AND STABILITY:

1. Store dried control material at 2-8°C (35-46°F). Stable to expiration date printed on the insert and label.
2. Store reconstituted control material at 2-8°C (35-46°F). See table for reconstituted stability.

V. PROCEDURE:

1. Remove cap from each vial to be used.
2. Reconstitute control material by adding exactly **10 mL of distilled water**, using a **10 mL** volumetric pipette or equivalent.
3. Replace cap and let sit 10-15 minutes.
4. Swirl gently 3-4 minutes to ensure a homogeneous mixture.
5. Swirl gently each time an aliquot is removed to ensure a homogeneous mixture.

VI. LIMITATIONS:

1. Results are dependent upon proper storage and adequate mixing.
2. Control material approximates a patient specimen.

VII. EXPECTED VALUES:

1. Laboratories should establish their own mean values; an individual laboratory's mean of several determinations may not duplicate the values listed below.

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UTAK Laboratories, Inc.