REF 50703

Product Description: DAU Low Cutoff 2

Product Number:

LOT D6410

Lot Number:

 2027-07-31

Expiration Date:

Analyte	Reconstituted Stability (Days)	Gravimetrically Weighed Value	Verified Value	Units
α-Hydroxyalprazolam	25	250	300	ng/mL
Amobarbital	25	250	246	ng/mL
d-Amphetamine	25	312.5	316	ng/mL
(±)-3,4-Methylenedioxyamphetamine (MDA)	25	312.5	318	ng/mL
d-Methamphetamine	25	312.5	323	ng/mL
(±)-3,4-Methylenedioxymethamphetamine (MDMA)	25	312.5	325	ng/mL
(±)-3,4-Methylenedioxyethylamphetamine (MDEA)	25	312.5	312	ng/mL
Benzoylcegonine	25	125	113	ng/mL
Butalbital	25	250	252	ng/mL
Codeine	25	375	449	ng/mL
Hydrocodone	25	125	118	ng/mL
Hydromorphone	25	125	135	ng/mL
Nordiazepam	25	250	310	ng/mL
Lorazepam	25	250	300	ng/mL
Methadone	25	250	290	ng/mL
Methadone Primary Metabolite (EDDP)	25	250	272	ng/mL
Methaqualone	25	250	236	ng/mL
Morphine-3-β-D-Glucuronide	25	608 (Equivalent to 375 free Morphine)	536.22 (Equivalent to 331 free Morphine)	ng/mL
6-Monoacetylmorphine (6-MAM)	25	12.5	14	ng/mL
Nortriptyline	25	250	207	ng/mL
Oxazepam	25	250	280	ng/mL
Oxycodone	25	62.5	64	ng/mL
Oxymorphone	25	62.5	68	ng/mL
Pentobarbital	25	250	242	ng/mL
Phencyclidine (PCP)	25	31.25	32	ng/mL
Phenobarbital	25	250	244	ng/mL
Norpropoxyphene	25	250	240	ng/mL
Secobarbital	25	250	238	ng/mL
Temazepam	25	250	309	ng/mL
(l)-9-Carboxy-11-nor-Delta-9-THC (THC-COOH)	25	18.75	18	ng/mL

Matrix: Urine
Product Form: Lyophilized
Storage: 2-8°C (35-46°F)
Vial Volume: 10 mL

ENGLISH:

I. INTENDED USE:

The 50703, DAU Low Cutoff 2, urine control is for use as a control material to evaluate the performance of an analytical test method.

II. SUMMARY AND PRINCIPLES:

Several different techniques are used for evaluating or estimating the variance of results. The three subjects summarized below must be considered with any test method.

1. PREVENTIVE MEASURES:

These measures are usually contained in the design of the test method and include consideration for reagents, equipment, and operator errors. These measures are designed to minimize variance.

2. CONTROL MEASURES:

When a control sample is analyzed at the same time and in the same manner as a patient specimen, an estimate of variance is obtained for the test method. This estimate of variance can be compared to the acceptable limits of variance of the test method.

3. STATISTICAL ANALYSIS OF PATIENT RESULTS:

As an aid in evaluating overall test results, the past experience of expected results can be compared to the results of any given test run. For example, it would not be expected that all results of a given test run be in an elevated range.

Control materials are widely used as a means to aid in the evaluation of test results. The following subjects are to be considered in the use of any control material.

1. Multi-Level	NORMAL / ELEVATED
2. Matrix	HUMAN / ANIMAL / CHEMICAL
3. Availability	SUFFICIENT FOR STATISTICS
4. Form	LIQUID / FROZEN / DRIED
5. Variety	DIFFERENT THAN CALIBRATORS

III. PRECAUTIONS:

1. Although the urine has been tested and found negative for HBsAg by RIA and HIV by EIA, the control material should be treated as any other potentially infectious agent.
2. For *in vitro* diagnostic use only.

IV. STORAGE AND STABILITY:

1. Store dried control material at 2-8°C (35-46°F). Stable to expiration date printed on the insert and label.
2. Store reconstituted control material at 2-8°C (35-46°F). See table for reconstituted stability.

V. PROCEDURE:

1. Remove cap from each vial to be used.
2. Reconstitute control material by adding exactly **10 mL of distilled water**, using a **10 mL** volumetric pipette or equivalent.
3. Replace cap and let sit 10-15 minutes.
4. Swirl gently 3-4 minutes to ensure a homogeneous mixture.
5. Swirl gently each time an aliquot is removed to ensure a homogeneous mixture.

VI. LIMITATIONS:

1. Results are dependent upon proper storage and adequate mixing.
2. Control material approximates a patient specimen.

VII. EXPECTED VALUES:

1. Laboratories should establish their own mean values; an individual laboratory's mean of several determinations may not duplicate the values listed below.

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