### ENGLISH:

#### I. INTENDED USE:

The Tricyclics VI MR, serum toxicology control is for use as a control material to evaluate the performance of an analytical test method.

#### II. SUMMARY AND PRINCIPLES:

Several different techniques are used for evaluating or estimating the variance of results. The three subjects summarized below must be considered with any test method.

##### 1. PREVENTIVE MEASURES:

These measures are usually contained in the design of the test method and include consideration for reagents, equipment, and operator errors. These measures are designed to minimize variance.

##### 2. CONTROL MEASURES:

When a control sample is analyzed at the same time and in the same manner as a patient specimen, an estimate of variance is obtained for the test method. This estimate of variance can be compared to the acceptable limits of variance of the test method.

##### 3. STATISTICAL ANALYSIS OF PATIENT RESULTS:

As an aid in evaluating overall test results, the past experience of expected results can be compared to the results of any given test run. For example, it would not be expected that all results of a given test run be in an elevated range.

Control materials are widely used as a means to aid in the evaluation of test results. The following subjects are to be considered in the use of any control material.

1. Multi-Level	NORMAL / ELEVATED
2. Matrix	HUMAN / ANIMAL / CHEMICAL
3. Availability	SUFFICIENT FOR STATISTICS
4. Form	LIQUID / FROZEN / DRIED
5. Variety	DIFFERENT THAN CALIBRATORS

#### III. PRECAUTIONS:

1. Although the serum has been tested and found negative for HBsAg by RIA and HIV by EIA, the control material should be treated as any other potentially infectious agent.
2. For *in vitro* diagnostic use only.

#### IV. STORAGE AND STABILITY:

1. Store dried control material at 2-8°C (35-46°F). Stable to expiration date printed on the insert and label.
2. Store reconstituted control material at or below 2-8°C (35-46°F). See table for reconstituted stability.

#### V. PROCEDURE:

1. Remove cap from each vial to be used.
2. Reconstitute control material by adding exactly **5 mL** of **distilled water**, using a **5 mL** volumetric pipette or equivalent.
3. Replace cap and let sit 10-15 minutes.
4. Swirl gently 3-4 minutes to ensure a homogeneous mixture.
5. Swirl gently each time an aliquot is removed to ensure a homogeneous mixture.

#### VI. LIMITATIONS:

1. Results are dependent upon proper storage and adequate mixing.
2. Control material approximates a patient specimen.

#### VII. EXPECTED VALUES:

1. Laboratories should establish their own mean values; an individual laboratory's mean of several determinations may not duplicate the values listed below.

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<b>REF</b>	66317	Product Description: Tricyclics VI MR		
Product Number:				
<b>LOT</b>	D3099	2025-04-30		
Lot Number: Expiration Date:				

Analyte	Stability(Days)	Gravimetrically Weighed Value	Verified Value	Units
Amitriptyline	25	100 <i>(0.36)</i>	80 <i>(0.28)</i>	ng/mL <i>(µmol/L)</i>
Doxepin	25	100 <i>(0.36)</i>	80 <i>(0.28)</i>	ng/mL <i>(µmol/L)</i>
N-Desmethyldoxepin	25	100 <i>(0.34)</i>	164 (64%) <i>(0.61)</i>	ng/mL <i>(µmol/L)</i>
Imipramine	25	100 <i>(0.36)</i>	86 <i>(0.30)</i>	ng/mL <i>(µmol/L)</i>
Desipramine	25	100 <i>(0.41)</i>	91 <i>(0.34)</i>	ng/mL <i>(µmol/L)</i>
Nortriptyline	25	100 <i>(0.38)</i>	104 <i>(0.39)</i>	ng/mL <i>(µmol/L)</i>

Matrix: Serum

Product Form: Lyophilized

Storage: 2-8°C (35-46°F)

Vial Volume: 5 mL

UTAK's express and implied warranties (including merchantability and fitness) are conditioned on the observance of UTAK's insert directions with respect to the use of UTAK's products.

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**For technical support please call: UTAK Technical Service (888) 882-5522**

**IVD**  
 UTAK Laboratories, Inc.  
 25020 Avenue Tibbitts  
 Valencia, CA 91355  
 T: (888) 882-5522  
 F: (661) 294-9272  
 E: [welovecontrol@utak.com](mailto:welovecontrol@utak.com)

**EMERGO EUROPE**  
 PRINSESSEGRACHT 20  
 2514 AP THE HAGUE  
 THE NETHERLANDS