REF 72741		Product Description: AED II MR		
Product Number:				
LOT D5368		2027-03-31		
Lot Number:		Expiration Date:		
Analyte	Reconstituted Stability (Days)	Gravimetrically Weighed Value	Verified Value	Units
Felbamate	25	50	49	µg/mL
Gabapentin	25	10	10.9	µg/mL
Lacosamide	25	5	4.8	µg/mL
Lamotrigine	25	5	3	µg/mL
Levetiracetam	25	15	15	µg/mL
Oxcarbazepine	25	2	Present ¹	µg/mL
10-OH-Carbamazepine	25	20	22	µg/mL
Pregabalin	25	5	4.9	µg/mL
Rufinamide	25	5	4.7	µg/mL
Tiagabine	25	80	85.5	ng/mL
Topiramate	25	15	13.6	µg/mL
Zonisamide	25	20	19	µg/mL
Matrix: Serum Product Form: Lyophilized Storage: 2-8°C (35-46°F) Vial Volume: 5 mL				

¹Analyte has been weighed and spiked into the control; however, no claim will be made for verified value.

ENGLISH:

I. INTENDED USE:

The 72741, AED II MR, serum control is for use as a control material to evaluate the performance of an analytical test method.

II. SUMMARY AND PRINCIPLES:

Several different techniques are used for evaluating or estimating the variance of results. The three subjects summarized below must be considered with any test method.

1. PREVENTIVE MEASURES:

These measures are usually contained in the design of the test method and include consideration for reagents, equipment, and operator errors. These measures are designed to minimize variance.

2. CONTROL MEASURES:

When a control sample is analyzed at the same time and in the same manner as a patient specimen, an estimate of variance is obtained for the test method. This estimate of variance can be compared to the acceptable limits of variance of the test method.

3. STATISTICAL ANALYSIS OF PATIENT RESULTS:

As an aid in evaluating overall test results, the past experience of expected results can be compared to the results of any given test run. For example, it would not be expected that all results of a given test run be in an elevated range.

Control materials are widely used as a means to aid in the evaluation of test results. The following subjects are to be considered in the use of any control material.

1. Multi-Level	NORMAL / ELEVATED
2. Matrix	HUMAN / ANIMAL / CHEMICAL
3. Availability	SUFFICIENT FOR STATISTICS
4. Form	LIQUID / FROZEN / DRIED
5. Variety	DIFFERENT THAN CALIBRATORS

III. PRECAUTIONS:

- Although the serum has been tested and found negative for HBsAg by RIA and HIV by EIA, the control material should be treated as any other potentially infectious agent.
- For *in vitro* diagnostic use only.

IV. STORAGE AND STABILITY:

- Store dried control material at 2-8°C (35-46°F). Stable to expiration date printed on the insert and label.
- Store reconstituted control material at 2-8°C (35-46°F). See table for reconstituted stability.

V. PROCEDURE:

- Remove cap from each vial to be used.
- Reconstitute control material by adding exactly **5 mL** of **distilled water**, using a **5 mL** volumetric pipette or equivalent.
- Replace cap and let sit 10-15 minutes.
- Swirl gently 3-4 minutes to ensure a homogeneous mixture.
- Swirl gently each time an aliquot is removed to ensure a homogeneous mixture.

VI. LIMITATIONS:

- Results are dependent upon proper storage and adequate mixing.
- Control material approximates a patient specimen.

VII. EXPECTED VALUES:

- Laboratories should establish their own mean values; an individual laboratory's mean of several determinations may not duplicate the values listed below.

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