REF 91002

Product Description: PM Plus High

Product Number:

LOT D6257

Lot Number:

 2027-06-30

Expiration Date:

Analyte	Reconstituted Stability (Days)	Gravimetrically Weighed Value	Verified Value	Units
Buprenorphine	25	400	268	ng/mL
Norbuprenorphine	25	400	256	ng/mL
Codeine	25	400	377	ng/mL
Norcodeine	25	400	273	ng/mL
Fentanyl	25	40	35	ng/mL
Norfentanyl	25	40	31	ng/mL
Hydrocodone	25	400	329	ng/mL
Norhydrocodone	25	400	285	ng/mL
Hydromorphone	25	400	312	ng/mL
Meperidine	25	400	435	ng/mL
Normeperidine	25	400	426	ng/mL
Methadone	25	400	323	ng/mL
Methadone Primary Metabolite (EDDP)	25	400	301	ng/mL
Morphine	25	400	319	ng/mL
6-Monoacetylmorphine (6-MAM)	25	40	40.2	ng/mL
Oxycodone	25	400	290	ng/mL
Noroxycodone	25	400	305	ng/mL
Oxymorphone	25	400	337	ng/mL
Noroxymorphone	25	400	228	ng/mL
Tapentadol	25	400	354	ng/mL
N-Desmethyltapentadol	25	400	403.1	ng/mL
Tramadol	25	400	341	ng/mL
O-Desmethyl-cis-Tramadol	25	400	329	ng/mL

Matrix: Urine
Product Form: Lyophilized
Storage: 2-8°C (35-46°F)
Vial Volume: 5 mL

ENGLISH:

I. INTENDED USE:

The 91002, PM Plus High, urine control is for use as a control material to evaluate the performance of an analytical test method.

II. SUMMARY AND PRINCIPLES:

Several different techniques are used for evaluating or estimating the variance of results. The three subjects summarized below must be considered with any test method.

1. PREVENTIVE MEASURES:

These measures are usually contained in the design of the test method and include consideration for reagents, equipment, and operator errors. These measures are designed to minimize variance.

2. CONTROL MEASURES:

When a control sample is analyzed at the same time and in the same manner as a patient specimen, an estimate of variance is obtained for the test method. This estimate of variance can be compared to the acceptable limits of variance of the test method.

3. STATISTICAL ANALYSIS OF PATIENT RESULTS:

As an aid in evaluating overall test results, the past experience of expected results can be compared to the results of any given test run. For example, it would not be expected that all results of a given test run be in an elevated range.

Control materials are widely used as a means to aid in the evaluation of test results. The following subjects are to be considered in the use of any control material.

1. Multi-Level	NORMAL / ELEVATED
2. Matrix	HUMAN / ANIMAL / CHEMICAL
3. Availability	SUFFICIENT FOR STATISTICS
4. Form	LIQUID / FROZEN / DRIED
5. Variety	DIFFERENT THAN CALIBRATORS

III. PRECAUTIONS:

1. Although the urine has been tested and found negative for HBsAg by RIA and HIV by EIA, the control material should be treated as any other potentially infectious agent.
2. For *in vitro* diagnostic use only.

IV. STORAGE AND STABILITY:

1. Store dried control material at 2-8°C (35-46°F). Stable to expiration date printed on the insert and label.
2. Store reconstituted control material at 2-8°C (35-46°F). See table for reconstituted stability.

V. PROCEDURE:

1. Remove cap from each vial to be used.
2. Reconstitute control material by adding exactly **5 mL** of **distilled water**, using a **5 mL** volumetric pipette or equivalent.
3. Replace cap and let sit 10-15 minutes.
4. Swirl gently 3-4 minutes to ensure a homogeneous mixture.
5. Swirl gently each time an aliquot is removed to ensure a homogeneous mixture.

VI. LIMITATIONS:

1. Results are dependent upon proper storage and adequate mixing.
2. Control material approximates a patient specimen.

VII. EXPECTED VALUES:

1. Laboratories should establish their own mean values; an individual laboratory's mean of several determinations may not duplicate the values listed below.

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