Reference Data Technical Document

This letter is to share further information on the 3rd party reference value that is indicated on the product insert for UTAK part number 50703 – DAU LOW CUTOFF 2 Urine Control, product lot number C6572.

This product was ran by an individual analysis by 3rd party laboratories that have unbiased methods and are ran as a patient samples to mimic a patient result. The results below are for reference only.

Analyte(s)	Gravimetric Spike Value	Verified Value	Method of Testing
Amobarbital	250 ng/mL	128 ng/mL	GC/MS
Lorazepam	250 ng/mL	313 ng/mL	LC/MS-MS
(I)-9-Carboxy-11-nor-Delta-9-THC (THC-COOH)	18.75 ng/mL	14 ng/mL	GC/MS

UTAK does not make claims regarding which lab results are correct or incorrect. This product is Class I unassayed quality control material. Each lab is to establish their own mean and acceptable range following their standard operating procedures related to the validation of quality control material.

Please contact our customer service team for further support or questions related to this product.

Approved by: Michelle Perdomo
QA/QC Specialist

Date: February 18, 2021

ENGLISH:

I. INTENDED USE:

The DAU Low Cutoff 2, urine toxicology control is for use as a control material to evaluate the performance of an analytical test method.

II. SUMMARY AND PRINCIPLES:

Several different techniques are used for evaluating or estimating the variance of results. The three subjects summarized below must be considered with any test method.

1. PREVENTIVE MEASURES:

These measures are usually contained in the design of the test method and include consideration for reagents, equipment, and operator errors. These measures are designed to minimize variance.

2. CONTROL MEASURES:

When a control sample is analyzed at the same time and in the same manner as a patient specimen, an estimate of variance is obtained for the test method. This estimate of variance can be compared to the acceptable limits of variance of the test method.

3. STATISTICAL ANALYSIS OF PATIENT RESULTS:

As an aid in evaluating overall test results, the past experience of expected results can be compared to the results of any given test run. For example, it would not be expected that all results of a given test run be in an elevated range.

Control materials are widely used as a means to aid in the evaluation of test results. The following subjects are to be considered in the use of any control material.

1. Multi-Level	NORMAL / ELEVATED
2. Matrix	HUMAN / ANIMAL / CHEMICAL
3. Availability	SUFFICIENT FOR STATISTICS
4. Form	LIQUID / FROZEN / DRIED
5. Variety	DIFFERENT THAN CALIBRATORS

III. PRECAUTIONS:

1. Although the urine has been tested and found negative for HBsAg by RIA and HIV by EIA, the control material should be treated as any other potentially infectious agent.
2. For *in vitro* diagnostic use only.

IV. STORAGE AND STABILITY:

1. Store dried control material at 2-8°C (35-46°F). Stable to expiration date printed on the insert and label.
2. Store reconstituted control material at or below 2-8°C (35-46°F). See table for reconstituted stability.

V. PROCEDURE:

1. Remove cap from each vial to be used.
2. Reconstitute control material by adding exactly **10 mL** of **distilled water**, using a 10 mL volumetric pipette or equivalent.
3. Replace cap and let sit 10-15 minutes.
4. Swirl gently 3-4 minutes to ensure a homogeneous mixture.
5. Swirl gently each time an aliquot is removed to ensure a homogeneous mixture.

VI. LIMITATIONS:

1. Results are dependent upon proper storage and adequate mixing.
2. Control material approximates a patient specimen.

VII. EXPECTED VALUES:

1. Laboratories should establish their own mean values; an individual laboratory's mean of several determinations may not duplicate the values listed below.

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UTAK Laboratories, Inc.