

AMISULPRIDE PLUS

SERUM TOXICOLOGY CONTROL

I. INTENDED USE:

Many analytes can be measured in patients by using analytical test methods. The UTAK Amisulpride Plus Control is for use as a quality control material for measuring the levels of analytes in serum. It is intended for use on a continuous basis so that a statistical evaluation of testing performance can be obtained.

II. SUMMARY AND PRINCIPLES:

Several different techniques are used for evaluating or estimating the variance of results. The three subjects summarized below must be considered with any test method.

1. PREVENTIVE MEASURES:

These measures are usually contained in the design of the test method and include consideration for reagents, equipment, and operator errors. These measures are designed to minimize variance.

2. QUALITY CONTROL MEASURES:

When a quality control sample is analyzed at the same time and in the same manner as a patient specimen, an estimate of variance is obtained for the test method. This estimate of variance can be compared to the acceptable limits of variance of the test method.

3. STATISTICAL ANALYSIS OF PATIENT RESULTS:

As an aid in evaluating overall test results, the past experience of expected results can be compared to the results of any given test run. For example, it would not be expected that all results of a given test run be in an elevated range.

Quality control materials are widely used as a means to aid in the evaluation of test results. The following subjects are to be considered in the use of any control material.

1. Multi-Level	NORMAL / ELEVATED
2. Matrix	HUMAN / ANIMAL / CHEMICAL
3. Availability	SUFFICIENT FOR STATISTICS
4. Form	LIQUID / FROZEN / DRIED
5. Variety	DIFFERENT THAN CALIBRATORS

The UTAK Amisulpride Plus Control is prepared from normal human materials and will generate data that checks and evaluates the results of a test method over the normal and elevated ranges. The principles of statistics require that the same material be available for comparison for any given time period. Dried control materials both extend the usable time period and allow larger quantities to be available. Statistical accuracy requires that a test method be defined for variance and be calibrated with a suitable standard. The quality control materials that are used must be of a sufficient variety so that the measurements and the data that are obtained are independent of the calibration standards. By using a variety of materials, the entire test method can be continuously evaluated to ensure reliable results.

III. PRODUCT DESCRIPTION:

The matrix for the UTAK Amisulpride Plus Control is prepared from normal human serum. Every analyte is added and adjusted to the desired concentration range for each lot prepared (Target Value). Quality control before, during, and after the preparation of the control material ensures that each lot is of the same quality.

IV. PRECAUTIONS:

- Although the serum has been tested and found negative for HBsAg by RIA and HIV by EIA, the control material should be treated as any other potentially infectious agent.
- For in vitro diagnostic use only.
- For analytical use only.

V. STORAGE AND STABILITY:

- Store dried control material at 2-8°C (35-46°F). Stable to expiration date printed on the insert and label.
- Store reconstituted control material at 2-8°C (35-46°F). **See table for reconstituted stability**

VI. PROCEDURE:

- Remove cap from each vial to be used.
- Reconstitute control material by adding **exactly 3.0 mL of distilled water**, using a **3.0 mL** volumetric pipette or equivalent.
- Replace cap and let sit 10-15 minutes.
- Swirl gently 3-4 minutes to ensure a homogeneous mixture.
- Swirl gently each time an aliquot is removed to ensure a homogeneous mixture.
- Assay control material in same manner as patient specimens, following the exact same instructions from the entire test method.
- Record the results obtained on a quality control chart that describes statistical limits for the test method and the particular lot of control material.

VII. LIMITATIONS:

- Control material is for use in quality control programs only; it is not intended for use as a calibration standard.
- Check the lot number on each vial to be sure it corresponds to the lot number printed on the insert.
- Results are dependent upon proper storage, reconstitution accuracy, and adequate mixing.
- Control material approximates a patient specimen.

VIII. EXPECTED VALUES:

- Laboratories should establish their own statistical values for precision and expected ranges; these values should fall within ±15% of the Target Value.

Lot Number : <u>A3366</u>		Expiration Date : <u>10/16</u>	
Analyte	Target Value	Units	Reconstituted Stability (days)
Amisulpride	200	µg/L	25
Amitriptyline	120	µg/L	30
Aripiprazole	200	µg/L	30
Atomoxetine	60	µg/L	25
Chlorprothixene	125	µg/L	25
Citalopram	100	µg/L	30
Clomipramine	120	µg/L	30
N-Desmethylclomipramine	120	µg/L	30
Clozapine	360	µg/L	30
N-Desmethylclozapine	180	µg/L	25
Desipramine	120	µg/L	25
Doxepin	100	µg/L	25
Desmethyldoxepin	100	µg/L	30
Duloxetine	20	µg/L	30
Fluoxetine	250	µg/L	30
Norfluoxetine	300	µg/L	25
Fluvoxamine	120	µg/L	25
Haloperidol	10	µg/L	30
Imipramine	120	µg/L	30
Levomepromazine	40	µg/L	30
Maprotiline	120	µg/L	25
Methylphenidate	20	µg/L	25
Mianserin	60	µg/L	30
Mirtazapine	150	µg/L	25
Nortriptyline	120	µg/L	25
Olanzapine	50	µg/L	30
N-Desmethylolanzapine	50	µg/L	30
Paroxetine	60	µg/L	30
Perazine	150	µg/L	25
Pipamperone	150	µg/L	30
Promethazine	150	µg/L	25
Protriptyline	120	µg/L	25
Quetiapine	180	µg/L	30
Reboxetine	40	µg/L	30
Risperidone	10	µg/L	25
9-Hydroxy Risperidone	50	µg/L	25
Ritalinic Acid	100	µg/L	30
Sertraline	60	µg/L	25
N-Desmethylsertraline	100	µg/L	30
Trazodone	1000	µg/L	30
Trimipramine	120	µg/L	25
N-Desmethyltrimipramine	120	µg/L	25
Venlafaxine	100	µg/L	25
O-Desmethylvenlafaxine	250	µg/L	30
Ziprasidone	50	µg/L	25

UTAK's express and implied warranties (including merchantability and fitness) are conditioned on the observance of UTAK's insert directions with respect to the use of UTAK's products.

For technical assistance call: UTAK Technical Service (800) 235-3442

UTAK LABORATORIES, INC.
 25020 AVENUE TIBBITTS
 VALENCIA, CA 91355
 TEL: (661) 294-3935
 FAX: (661) 294-9272
 E-MAIL: INQUIRIES@UTAK.COM

PRODUCT NUMBER:

 1x3ML VIALS, DRIED

EC AUTHORIZED REPRESENTATIVE
 EMERGO EUROPE
 MOLENSTRAAT 15
 2513 BH, THE HAGUE
 THE NETHERLANDS

