UTAK

DAU High Cutoff 1
Urine Control

REF 50700
Product Number:

Product Description: DAU High Cutoff 1

LOT D9113
Lot Number:

2029-03-31
Expiration Date:

Analyte	Reconstituted Stability (Days)	Gravimetrically Weighed Value	Verified Value	Units
α-Hydroxyalprazolam	25	225	233	ng/mL
Amobarbital	25	225	215	ng/mL
d-Amphetamine	25	375	382	ng/mL
(±)-3,4-Methylenedioxyamphetamine (MDA)	25	375	368	ng/mL
d-Methamphetamine	25	375	360	ng/mL
(±)-3,4-Methylenedioxymethamphetamine (MDMA)	25	375	372	ng/mL
(±)-3,4-Methylenedioxyethylamphetamine (MDEA)	25	375	405	ng/mL
Benzoylecgonine	25	112.5	110	ng/mL
Butalbital	25	225	203	ng/mL
Codeine	25	1500	1469	ng/mL
Hydrocodone	25	300	266	ng/mL
Hydromorphone	25	300	295	ng/mL
Nordiazepam	25	225	234	ng/mL
Lorazepam	25	225	225	ng/mL
Methadone	25	225	205	ng/mL
Methadone Primary Metabolite (EDDP)	25	225	201	ng/mL
Methaqualone	25	225	228	ng/mL
Morphine-3-β-D-Glucuronide	25	2430 (Equivalent to 1500 free Morphine)	2304 (Equivalent to 1422 free Morphine)	ng/mL
6-Monoacetylmorphine (6-MAM)	25	7.5	7	ng/mL
Nortriptyline	25	225	181	ng/mL
Oxazepam	25	225	231	ng/mL
Oxycodone	25	225	217	ng/mL
Oxymorphone	25	225	240	ng/mL
Pentobarbital	25	225	204	ng/mL
Phencyclidine (PCP)	25	18.75	19	ng/mL
Phenobarbital	25	225	222	ng/mL
Norpropoxyphene	25	225	178	ng/mL
Secobarbital	25	225	219	ng/mL
Temazepam	25	225	224	ng/mL
(I)-9-Carboxy-11-nor-Delta-9-THC (THC-COOH)	25	11.25	11	ng/mL

Matrix: Urine
Product Form: Lyophilized
Storage: 2-8°C (35-46°F)
Vial Volume: 10 mL

ENGLISH:

I. INTENDED USE:
The 50700, DAU High Cutoff 1, urine control is for use as a control material to evaluate the performance of an analytical test method.

II. SUMMARY AND PRINCIPLES:
Several different techniques are used for evaluating or estimating the variance of results. The three subjects summarized below must be considered with any test method.

1. PREVENTIVE MEASURES:
 These measures are usually contained in the design of the test method and include consideration for reagents, equipment, and operator errors. These measures are designed to minimize variance.
2. CONTROL MEASURES:
 When a control sample is analyzed at the same time and in the same manner as a patient specimen, an estimate of variance is obtained for the test method. This estimate of variance can be compared to the acceptable limits of variance of the test method.
3. STATISTICAL ANALYSIS OF PATIENT RESULTS:
 As an aid in evaluating overall test results, the past experience of expected results can be compared to the results of any given test run. For example, it would not be expected that all results of a given test run be in an elevated range.

Control materials are widely used as a means to aid in the evaluation of test results. The following subjects are to be considered in the use of any control material.

1. Multi-Level	NORMAL / ELEVATED
2. Matrix	HUMAN / ANIMAL / CHEMICAL
3. Availability	SUFFICIENT FOR STATISTICS
4. Form	LIQUID / FROZEN / DRIED
5. Variety	DIFFERENT THAN CALIBRATORS

III. PRECAUTIONS:

1. Although the urine has been tested and found negative for HBsAg by RIA and HIV by EIA, the control material should be treated as any other potentially infectious agent.
2. For *in vitro* diagnostic use only.

IV. STORAGE AND STABILITY:

1. Store dried control material at 2-8°C (35-46°F). Stable to expiration date printed on the insert and label.
2. Store reconstituted control material at 2-8°C (35-46°F). See table for reconstituted stability.

V. PROCEDURE:

1. Remove cap from each vial to be used.
2. Reconstitute control material by adding exactly **10 mL** of **distilled water**, using a **10 mL** volumetric pipette or equivalent.
3. Replace cap and let sit 10-15 minutes.
4. Swirl gently 3-4 minutes to ensure a homogeneous mixture.
5. Swirl gently each time an aliquot is removed to ensure a homogeneous mixture.

VI. LIMITATIONS:

1. Results are dependent upon proper storage and adequate mixing.
2. Control material approximates a patient specimen.

VII. EXPECTED VALUES:

1. Laboratories should establish their own mean values; an individual laboratory's mean of several determinations may not duplicate the values listed below.

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UTAK Laboratories, Inc.