



people are dying.
your results save lives.

fentanyl analogues

quality controls in whole blood & urine



UTAK

fentanyl analogues

4-ANPP	m-Methylmethoxyacetylfentanyl
Acetylfentanyl	Norfentanyl
Acrylfentanyl	Ocfentanyl
Alfentanil	ortho-Fluorofentanyl
Butyryl Fentanyl	para-Fluorobutyrylfentanyl
Carfentanyl	para-Fluorofentanyl
cis-3-Methylfentanyl	para-Fluoroisobutyrylfentanyl
Crotonylfentanyl	Sufentanil
Cyclopropylylfentanyl	Tetrahydrofuran Fentanyl
Fentanyl	U-47700
Furanyl Fentanyl	U-49900
Isobutyrylfentanyl	U51754
Methoxyacetylfentanyl	Valerylfentanyl

matrix	whole blood	urine
concentration	2 ng/mL	5 ng/mL
unit size	5x3mL	5x3mL
form	frozen	frozen
shelf life	18 months	18 months
stability	25 days	25 days

Custom concentrations also available.

Get a sample: utak.com/fentanyl

UTAK



25020 Avenue Tibbitts
Valencia, CA 91355
utak.com

T 661.294.3935
F 661.294.9272