



November 11, 2019

UTAK Product Change Notification

UTAK Laboratories, Inc has made a design change to UTAK part number 17013 – Validity level 4. The design change of part number 17013 – Validity level 4 is specific to the pH concentration in the product. The design change to this product is effective with the newest manufactured lot # C4106. Please see the specific design change below.

The pH for the part #17013 has been changed from 4.5 – 5 to 2.0 – 2.8

For additional information, support or concerns related to this technical bulletin, please do not hesitate to contact UTAK Laboratories, Inc at 800-235-3442.

Sincerely,

Mathew Kopp
Chief Marketing Officer
UTAK Laboratories, Inc

